

## TruBio™ Product Family

Features	TruBio™ OS	TruBio™ SCADA	TruBio™ Monitor
<b>Basic modules and graphics</b>	Yes	SCADA loop blocks	Read only, no control
<b>TruConfig™</b> <ul style="list-style-type: none"> <li>• Config Builder</li> <li>• Save/Load</li> <li>• Signal Conditioning</li> </ul>	Yes	Yes	No Freeform Calc No Dropdown list
<b>TruConnect™</b> <ul style="list-style-type: none"> <li>• Functions can be restricted according to the TruBio licenses</li> </ul>	Yes	Yes	Read field values only No control
<b>TruData™</b> (15 points/vessel default) <ul style="list-style-type: none"> <li>• Offline Data Entry</li> </ul>	Data entry with auto standardize	Data entry only	Data entry only
<b>TruCal™</b> <ul style="list-style-type: none"> <li>• Text Prompts Guide</li> <li>• Electronic Report</li> </ul>	Yes	N/A	N/A
<b>TruHealth™</b> <ul style="list-style-type: none"> <li>• Secondary Measurement Signals</li> <li>• Predictive Measurement Control</li> </ul>	Yes	If data is available	If data is available
<b>Simulate</b>	Yes	N/A	N/A
<b>Snap-On</b> <ul style="list-style-type: none"> <li>• CIP/SIP skid</li> <li>• Batch Lite</li> <li>• XLReporter</li> <li>• PI</li> </ul>	Able to apply	Able to apply	No Batch Lite

