

TruBatch™

The TruBio™ OS has been designed to provide seamless integration with the S88 compliant suite of DeltaV and OSIsoft Batch software. The suite of DeltaV and OSIsoft Batch software products includes the following:

Batch Software Products

Basic Batch (Emerson/Delta V)	Entry level Batch product consists of the DeltaV Batch Executive (upgradeable)
BatchLite (Finesse)	Entry level Batch solution from TruBio™ SnapOn option
Advanced Batch (Emerson Delta V)	Mid-level Batch product consists of the DeltaV Batch Executive, Batch Historian and Soft Phase Server (upgradable)
Professional Batch (Emerson Delta V)	High-level Batch product consists of the DeltaV Batch Executive, Batch Historian, Soft Phase Server, Campaign Manager and Recipe Exchange
PI Batch (OSIsoftPI)	Highly capable Batch data management tool can work alongside the DeltaV Batch suite of products

It is important to note that the use of either the DeltaV or OSIsoft PI Batch software requires the purchase of additional licensing and custom configuration. For assistance with specifying Batch licensing or to discuss your custom configuration requirements, please contact your local Finesse representative.



Features

- Coordinates recipe execution
- Performs equipment arbitration
- Archives recipe execution and event data
- Supports batch based analysis
- Able to run dozens of recipes simultaneously and several hundred phases