



## TruCables™

The TruCables™ family of products is designed for improved performance, better repeatability and increased reliability in bioprocess applications such as cell culture and fermentation.

### Specifications

pH	S8, K8 and VP connectors
DO	D4, T82 and VP connectors
Lengths	6 ft (2 m) and 10 ft (3 m)

### pH Cables

S8	K8	VP
CAB-S8-06, 6 ft/2 m	CAB-K8-06, 6 ft/2 m	CAB-VP-06, 6 ft/2 m
CAB-S8-10, 10 ft/3 m	CAB-K8-10, 10 ft/3 m	CAB-VP-10, 10 ft/3 m

### DO Cables

D4/T82	VP
CAB-D4-06, 6 ft/2 m	CAB-VP-06, 6 ft/2 m
CAB-D4-10, 10 ft/3 m	CAB-VP-10, 10 ft/3 m

## Features


- Optimum performance
- Designed for high performance
- Compatible with Broadley-James and Mettler Toledo
- Available with S8, K8, D4/T82 or VP connectors
- Lengths of 6ft (2 m) or 10 ft (3 m)




S8



K8



D4/T82



VP