

Does your Bioreactor Controller have the ability to monitor cell growth real time – and change feed strategies accordingly?

Now it can...

Finesse Solutions

Introducing *TruCell2* 
the next generation in on-line Cell Density measurement

TruCell2™ is the life science industry's first and only cell density solution that is both easy to use and has a cost comparable to Dissolved Oxygen. The TruCell2™ solution from Finesse delivers increased performance, repeatability and reliability in a user friendly design.

Best in Class performance means the TruCell2™ is able to deliver optical measurements with unequalled Precision, Zero Drift and Dynamic Range.

Intelligent Transmitter design allows you to easily enable Automatic Zeroing, Curve Fitting, and for the first time, Temperature Compensation.

Improved Sensor design delivers 50+ auto-claves @ 134°C, USP Class VI Wetted Materials and Optical Filters for sensor verification.

Contact Us To learn more about Finesse and the TruCell2™ give us a call or visit us on the web. Our team of industry professionals is here to assist you.

